

# Chemistry

## Laws of gases



### ■ C11 - Determination of n: quantity of matter in a gas

#### ■ Theoretical notes and objectives

The law of perfect gases shows that the pressure P of a gas, the volume V, the quantity of matter n and the absolute temperature T are related by the equation:


$$P V = n R T$$

Where  $n = \frac{PV}{RT}$

The relation between the temperature t expressed in °C and the absolute temperature T is written as:

$$T = 273.15 + t$$

The reaction studied here is written as:



An excess of nitric acid shall be used to make sure that the reaction is completed.

Aim of the experiment:

Study the law linking the pressure P and the volume V, the quantity of matter n and the absolute temperature T.


Obtain an estimate of the quantity of matter n.

#### ■ Experimental procedure

Weight about 2g of sodium bicarbonate and place it in the flask.

Prepare 10 mL of nitric acid at 6 mol/L and place it in a small beaker.

#### ■ Setting the apparatus







#### ■ Material


VTT console	ref 477000
VTT pressure sensor	ref 472010
VTT thermometer sensor	ref 472004
500mL glass flask	ref 713621
Stopper for the flask	ref 253082
Beaker (25 mL) - pack of 10	ref 713502
Pipette (10 mL)	ref 713047
Dropper bulb	ref 703172
Round wash bottles (pack of 6)	ref 723208
Distilled water (1L)	ref 107340
Nitric acid 62% *	ref 106012
Sodium bicarbonate 250g	ref 107270
Balance	ref:701276
Generis 5+ software	ref 000720

\* Be careful, dangerous product. Handle with care.




## ■ Settings required before the experiment

Work to be done	How to do it?
<p>Connect the VTT interface to the computer. Calibration of sensors.</p> <p>Launch Generis</p> <p>Put the adapters in ordinate. Give a name to the magnitude.</p> <p>Pressure sensor in channel /.</p> <p>Temperature sensor in channel //.</p> <p>Time in abscissa. Duration: 180 s</p> <p>Start the acquisition: Tilt the flask pouring the nitric acid on the bicarbonate. The reaction is over when there are no more gas bubbles. Observe the increase in temperature and pressure</p> <p>Save the file corresponding to this graph <math>P = f(V)</math> under the name <b>Determination of n.lab</b>.</p>	<p>Refer to the sensors' instruction manuals. Start, Programmes, Science Studio, Generis 5+</p> <p>Drag-drop the <b>pressure sensor</b>  on one of the channels in ordinate and name the magnitude P.</p> <p>Drag-drop the <b>temperature sensor</b>  on one of the channels in ordinate and name the magnitude <math>\theta</math>. (use the F12 key for the Greek alphabet). Unit: °C</p> <p>Drag-drop the <b>icon time</b>  in abscissa. Name the magnitude t. Duration: 180 s</p> <p>Click on </p> <p>Go to File, Save as... Select the directory folder and give the name <b>Determination of n</b>.</p>

## ■ Results obtained



## ■ Analysis of results

Work to be done	How to do it?
<p>Magnify the curves: use the automatic scales</p> <p>Using the relative pointer, determine the increase in pressure and the rise in temperature.</p> <p>Save the work completed</p>	<p>Click on the icon </p> <p>Right click on the graph and select <b>Pointer</b>. Click on the first summit and without dropping go to the next summit. The value of increase in pressure can be read in the caption bar.</p> <p>Change the axis for reading the second value: icon </p> <p>Click on the icon  or Menu File, Save.</p>

## ■ Interpretation of results

Note the increase of pressure and the rise in temperature.






Explain the increase of pressure.

Calculate the number of moles of bicarbonate used knowing the pressure of CO<sub>2</sub> gas, its temperature and the volume of calorimeter: 500mL

Calculate the number of moles of bicarbonate used based on the weight indicated when starting (4g).

Compare the two results. Explain the difference between the two results.

## ■ Drafting of report


Work to be done	How to do it
<p>Click on the tab Report</p> <p>Enter the title</p> <p>Insert the graph P(t) and <math>\theta(t)</math>: a single graph</p> <p>Answer the questions asked.</p> <p>Go to Print preview and if required, modify the elements of the report.</p>	<p></p> <p>Click on the icon  and define a rectangle Enter the title</p> <p>Open the dialog box  and click on the curves to be inserted P(t) and <math>\theta(t)</math>. Define a rectangle containing the curves.</p> <p>Click on the icon  and define a rectangle for the answers.</p> <p>Click on the icon </p>

# Chemistry

## Laws of gases



### ■ Determination of n: quantity of matter in a gas – Correction sheet



Increase of the pressure:  $\Delta P = 816 \text{ hPa} = 81600 \text{ Pa}$

Rise in the temperature:  $\Delta t = 0.75^\circ\text{C}$ . Final temperature  $t = 24^\circ\text{C}$ .

The increase of pressure is linked to the release of gaseous carbon dioxide in the flask.

Number of moles of bicarbonate used:

$$n = \frac{PV}{RT} = \frac{82000 \cdot 0,0005}{8.32 \cdot 297} = 0.017 = 1.7 \cdot 10^{-2} \text{ mol}$$

Number of moles of  $\text{NaHCO}_3$  used (weighing):

$$n = 1.90 / 84.01 = 2.26 \cdot 10^{-2} \text{ mol}$$

For obtaining more precise results, a sealed calorimeter should be used.